




REHASENSE



HEAVY DUTY



SERVER HD

THE HEAVY DUTY ROLLATOR

The Heavy Duty version of the Rehasense Server rollator range. Weighing less than 8 kg, the Server HD re-enforced, extra wide frame can carry 200 kg.

SERVER HD

THE HEAVY DUTY ROLLATOR



REHASENSE

FOLDABLE AND EASY TO CARRY

Server HD folds easily and stays that way with a simple locking system that doubles up as an ergonomically shaped carrying handle.

STEADY AND DURABLE

Frame and seat are tested and approved for maximum user weight of 200 kg. The brake mechanism is light yet positive. TPE-covered, unbreakable wheels.

HEIGHT ADJUSTABLE PUSH HANDLES

Server HD push handles are height adjustable with 10 positions, from 74 cm to 102 cm.


The seat height is 62 cm with an extra wide 55 cm seat base. User comfort is ensured with an optimal frame design and ergonomic push handles.

ACCESSORIES

Server HD can be fitted with a broad range of Rehasense accessories.



Smooth adjustable handle height and lockable handbrakes.



Seat height 62 cm.



Maintenance free brakes and TPE-covered unbreakable wheels.



Unique reinforced triangle profile and extra cross bar re-enforcement ensures to carry weight till 200 kg.



Stays firmly closed with the lock and is easy to open with the push button.



Practical easy-to-remove shopping bag with special designed easy-to-go clips.

TECHNICAL INFORMATION

Max. user weight: 200 kg

Item no:

Large: HRSGL600

Medium: HRSGM550

Measurements:

(width x depth x height)

Large: 69 x 68 x 74-102 cm

Medium: 69 x 68 x 66 - 86 cm

Seat width: 55 cm

Transportation measurements folded:

(depth x height x width)

Large: 68 x 86 x 22,5 cm

Medium: 68 x 79 x 22,5 cm

Weight:

Large: 7,9 kg

Medium: 7,6 kg

Shopping bag:

(width x depth x height)

35 x 17 x 25 cm

Max. load 5 kg.

Colour:

Sober Grey

Rehasense® Europe

Rehasense Sp. z o.o.

Sulejowska 45

97-300 Piotrków Trybunalski

Poland

Rehasense Benelux

Markeweg 81

8398 GM Blesdijke

Netherlands

Rehasense Deutschland

Briener Strasse 25

47533 Kleve

Deutschland



Extra wide seat 55 cm.

SERVER HD is CE-marked and fulfills the requirements of EU directive 93/42 EEC for medical devices.

Copyright © Rehasense® 2015