

SERVER

THE ELEGANT CLASSIC



REHASENSE



SERVER

THE ELEGANT CLASSIC



THE BESTSELLER

The SERVER is the classic aluminum lightweight rollator. The timelessly elegant frame is inspired by the purist style of Danish Design. Thanks to the unique profile the rollator can be easily kept clean and maintained. The ergonomic handle grips under the „Click & Safe“ safety catch allow secure lifting, carrying and transporting of the rollator. In the XL version the SERVER is available as the SERVER HD.

FEATURES AT A GLANCE

- Large shopping bag for carrying up to 5 kg
- User friendly and easily adjustable handle height
- Unique carrying handle: ergonomic handle grips under the mounting bars for lifting and carrying the rollator
- All-round safety reflectors available for shopping bag, frame and handles
- Large selection of accessories such as backrest straps, bags, tray, one-handed brake operation, slowdown brake, rollator umbrella

SIZES & COLOURS


Large: Champagne, Wine red, Brown, Anthracite, Black, White
Medium: Champagne, Wine red, Anthracite
Small: Champagne



TECHNICAL DETAILS

Weight	Max. load capacity	Seat height	Seat width	Width folded	Total width	Distance between handles	Handle height from small to large:
Large 7.0 kg	up to 150 kg	Large 62 cm	46 cm	23 cm	61 cm	47 cm	63-102 cm
Medium 6.8 kg		Medium 55 cm					
Small 6.7 kg excl. shopping bag 340 grams		Small 50 cm					

SPECIAL FEATURES



With classic TPE or shock-absorbent soft wheels

Including cane holder

Ergonomic grips moulded to shape

Brake lever can be operated without much effort

Ergonomic handle grips for lifting and carrying

Proven and secure "Click & Safe" locking system


ACCESSORIES



Back strap

LED-Lamps

Comfort seat



Rehasense's quality management complies with both the ISO 9001 and ISO 13485 standards, as confirmed by TÜV NORD. Rehasense's medical devices are CE-marked and meet the requirements of EU Regulation 2017/745 (MDR).

